

# The Oncomine™ BRCA Expanded next-generation sequencing assay: Development and analytical validation of a new panel for detection of SNV, insertions, deletions and copy number variants in a panel of 15 genes involved in homologous recombination repair of double-strand break DNA damage

Charles Scafe, Chenchen Yang, Yun Zhu, Yu-ting Tseng, Brooke McKnight, Fernando Farfan, Santhoshi Bandla, Seth Sadis, Steve Roman, Fiona C.L. Hyland  
Thermo Fisher Scientific, 200 Oyster Point Boulevard, South San Francisco, CA; Carlsbad, CA; Ann Arbor, MI

## INTRODUCTION

Breast cancer is the second most commonly diagnosed cancer (11.6% of total cases) world wide<sup>1</sup>. Many tumors display phenotypic similarity across tissue types, characterized by defects in biological pathways leading to susceptibility to particular therapies<sup>2</sup>. Deficiencies in one such pathway, homologous recombination DNA damage repair (HR DDR), affect the response of breast and ovarian cancers to agents that inhibit PARP function.

The Oncomine™ BRCA Expanded panel targets 15 genes (BRCA1 and 2, ATM, BARD1, BRIP1, CDK12, CHEK2, FANCD2, MRE11, NBN, PALB2, PPP2R2A, RAD51B, RAD54L, and TP53) that play various roles in the HR DDR pathway. The assay includes 1011 amplicons covering both hotspot variants and coding regions of these genes. An additional set of 255 pre-qualified genes is available for assay customization, including a recommended set of 25 genes of particular relevance to HR DDR pathway.

Small variants (SNV and INDEL) and large rearrangements (LRs) such as copy number variations (CNV) in these genes can be detected with next-generation sequencing (NGS). This panel is designed for FFPE samples and can also be used with samples having non-degraded DNA. Germline and somatic mutations, such as SNPs, Indels, Copy number gains, can be detected for all genes, and exon deletions can be detected for *BRCA1* and *BRCA2*.

## MATERIALS AND METHODS

The panel uses Ion AmpliSeq technology and is run on the Ion S5 Gene Studio™ instrument. Automated variant calling workflow is provided within Ion Reporter™. AcroMetrix Oncology Hotspot Control (AOHC) samples are used for verification of small variants calling. Samples with known gene copy number gains are used for CNV calling verification.

## RESULTS

### Genes in Oncomine™ BRCA Expanded panel


Table 1. 15 gene panel

| Gene    | Num_Amps | Coverage |
|---------|----------|----------|
| ATM     | 194      | 99.6     |
| BARD1   | 42       | 100      |
| BRCA1   | 113      | 100      |
| BRCA2   | 152      | 100      |
| BRIP1   | 72       | 99.9     |
| CDK12   | 67       | 100      |
| CHEK2   | 29       | 99.3     |
| FANCD2  | 81       | 98.4     |
| MRE11   | 42       | 100      |
| NBN     | 38       | 100      |
| PALB2   | 60       | 100      |
| PPP2R2A | 26       | 100      |
| RAD51B  | 23       | 100      |
| RAD54L  | 39       | 100      |
| TP53    | 23       | 100      |

Table 2. 25 expanded genes

| Gene     | Num_Amps | Coverage |
|----------|----------|----------|
| CHEK1    | 26       | 100      |
| EPCAM    | 27       | 96.6     |
| FANCB    | 35       | 100      |
| FANCC    | 34       | 100      |
| FANCE    | 26       | 99.2     |
| FANCF    | 11       | 100      |
| FANCG    | 30       | 100      |
| FANCI    | 72       | 100      |
| FANCL    | 27       | 100      |
| FANCM    | 84       | 100      |
| GEN1     | 37       | 99.4     |
| MLH3     | 56       | 100      |
| PARP1    | 54       | 100      |
| PMS1     | 37       | 98.7     |
| PRKDC    | 204      | 99.6     |
| RAD50    | 63       | 99.5     |
| RAD51    | 20       | 93.8     |
| RAD51C   | 23       | 100      |
| RAD51D   | 20       | 100      |
| RAD52    | 26       | 100      |
| RNASEH2A | 14       | 100      |
| RNASEH2B | 22       | 100      |
| RNASEH2C | 10       | 98.9     |
| TP53BP1  | 84       | 99.9     |
| XRCC2    | 11       | 100      |

Figure 1. Customized panel can include any of 270 genes including expanded genes, and can remove genes in the core panel



## General performance evaluation

Figure 2. Amplicon design and sequencing results






Figure 3. On target and uniformity



Average percent reads on target and uniformity of amplicon coverage of Oncomine™ BRCA Expanded panel on non-FFPE, good-quality FFPE and poor-quality FFPE samples.


Figure 4. Sensitivity and specificity in variants detection



Sensitivity and specificity of Oncomine™ BRCA Expanded panel based on 144 variants detection results using AcroMetrix Oncology Hotspot Control (AOHC) samples.

## Panel Selection on AmpliSeq.com

Figure 5. Oncomine™ tumor specific panels on AmpliSeq.com



## Variant calling examples

Figure 6. A gene with known copy number gains is being detected





Table 3. 32 out of 33 SNVs in TP53 are being detected using AOHC samples

| Locus         | COSMIC ID   | Genotype | Ref | Function   | Protein Change | Coverage | Alt Allele Frequency | Detected |
|---------------|-------------|----------|-----|------------|----------------|----------|----------------------|----------|
| chr17:7572986 | COSM307348  | G/A      | G   | nonsense   | p.Gln375Ter    | 1498     | 9.95                 | Yes      |
| chr17:7573010 | COSM1191161 | T/C      | T   | unknown    | NA             | 1492     | 10.19                | Yes      |
| chr17:7574003 | COSM11073   | G/A      | G   | nonsense   | p.Arg342Ter    | 2000     | 3.2                  | No       |
| chr17:7574012 | COSM11286   | C/A      | C   | nonsense   | p.Glu339Ter    | 1999     | 6.25                 | Yes      |
| chr17:7574018 | COSM11071   | G/A      | G   | missense   | p.Arg337Cys    | 1999     | 6.35                 | Yes      |
| chr17:7574026 | COSM11514   | C/A      | C   | missense   | p.Gly334Val    | 1998     | 6.31                 | Yes      |
| chr17:7576855 | COSM11354   | G/A      | G   | nonsense   | p.Gln331Ter    | 1997     | 6.31                 | Yes      |
| chr17:7576865 | COSM44823   | A/C      | A   | nonsense   | p.Tyr327Ter    | 2000     | 6.5                  | Yes      |
| chr17:7576883 | COSM46088   | T/C      | T   | synonymous | p.Lys521=      | 2000     | 6.4                  | Yes      |
| chr17:7576897 | COSM10786   | G/A      | G   | nonsense   | p.Gln317Ter    | 1996     | 6.71                 | Yes      |
| chr17:7577022 | COSM10663   | G/A      | G   | nonsense   | p.Arg306Ter    | 1717     | 6.99                 | Yes      |
| chr17:7577046 | COSM10710   | C/A      | C   | nonsense   | p.Glu298Ter    | 1720     | 6.86                 | Yes      |
| chr17:7577105 | COSM10863   | G/A      | G   | missense   | p.Pro278Leu    | 2000     | 6.55                 | Yes      |
| chr17:7577120 | COSM10660   | C/T      | C   | missense   | p.Arg273His    | 2000     | 5.95                 | Yes      |
| chr17:7577338 | COSM10662   | C/T      | C   | missense   | p.Arg248Gln    | 2000     | 5.25                 | Yes      |
| chr17:7577548 | COSM6932    | C/T      | C   | missense   | p.Gly245Ser    | 2000     | 7.4                  | Yes      |
| chr17:7577559 | COSM10812   | G/A      | G   | missense   | p.Ser241Phe    | 2000     | 6.9                  | Yes      |
| chr17:7577580 | COSM10725   | T/C      | T   | missense   | p.Tyr234Cys    | 2000     | 6.6                  | Yes      |
| chr17:7578190 | COSM10758   | T/C      | T   | missense   | p.Tyr220Cys    | 2000     | 6.2                  | Yes      |
| chr17:7578196 | COSM44317   | A/T      | A   | missense   | p.Val218Glu    | 1997     | 6.11                 | Yes      |
| chr17:7578203 | COSM10667   | C/T      | C   | missense   | p.Val216Met    | 2000     | 6.25                 | Yes      |
| chr17:7578235 | COSM43947   | T/C      | T   | missense   | p.Tyr205Cys    | 1998     | 6.06                 | Yes      |
| chr17:7578388 | COSM10738   | C/T      | T   | missense   | p.Arg181His    | 1942     | 7.21                 | Yes      |
| chr17:7579295 | COSM44985   | C/T      | C   | unknown    | NA             | 2000     | 5.65                 | Yes      |
| chr17:7579312 | COSM43904   | C/T      | C   | synonymous | p.Thr125=      | 2000     | 5.4                  | Yes      |
| chr17:7579358 | COSM10716   | C/A      | C   | missense   | p.Arg110Leu    | 2000     | 6.05                 | Yes      |
| chr17:7579368 | COSM46103   | A/C      | A   | missense   | p.Tyr107Asp    | 1999     | 6.15                 | Yes      |
| chr17:7579472 | COSM250611  | G/C      | G   | missense   | p.Pro72Arg     | 1999     | 94.2                 | Yes      |
| chr17:7579521 | COSM132168  | C/A      | C   | missense   | p.Glu56Ter     | 2000     | 5.95                 | Yes      |
| chr17:7579536 | COSM44907   | C/A      | C   | nonsense   | p.Glu51Ter     | 1998     | 6.01                 | Yes      |
| chr17:7579553 | COSM43664   | A/G      | A   | missense   | p.Leu45Pro     | 1999     | 5.7                  | Yes      |
| chr17:7579575 | COSM46286   | G/A      | G   | nonsense   | p.Gln38Ter     | 1996     | 5.96                 | Yes      |
| chr17:7579801 | NA          | G/C      | G   | unknown    | NA             | 1991     | 88.5                 | Yes      |

Figure 7. Oncomine™ tumor specific panel exon deletions detected in germline and FFPE samples



## Additional Oncomine™ tumor-specific panels

Figure 8. List of additional Additional Oncomine™ tumor-specific panels

Ion Torrent™ Oncomine™ Bladder Panel

Ion Torrent™ Oncomine™ Kidney Panel

Ion Torrent™ Oncomine™ BRCA Expanded Panel

Ion Torrent™ Oncomine™ Liver Panel

Ion Torrent™ Oncomine™ Colorectal and Pancreatic Panel

Ion Torrent™ Oncomine™ Lymphoma Panel

Ion Torrent™ Oncomine™ Gastric and Esophageal Panel

Ion Torrent™ Oncomine™ Melanoma Panel

Ion Torrent™ Oncomine™ Gynecological Panel

Ion Torrent™ Oncomine™ Prostate Panel

## HIGHLIGHTS

- Customizable panel with verified performance for clinical research
- Low sample input requirement (20 ng input DNA)
- Works with FFPE tissue samples
- End-to-end workflow including bioinformatics and reporting solutions
- Quick turnaround time (2 days from DNA to data)

## CONCLUSIONS

An NGS assay with a comprehensive data analysis approach was developed that is capable of detecting both small mutations and LR simultaneously in FFPE samples with high sensitivity. This is an important assay for *BRCA1/2* and HR DDR translational research.

## REFERENCES

- Bray, F. , Ferlay, J. , Soerjomataram, I. , Siegel, R. L., Torre, L. A. and Jemal, A. (2018), Global cancer statistics 2018: GLOBOCAN estimates of incidence and mortality worldwide for 36 cancers in 185 countries. CA: A Cancer Journal for Clinicians, 68: 394-424. doi:10.3322/caac.21492
- Heeke, A. L., Pishvaian, M. J., Lynce, F., Xiu, J., Brody, J. R., Chen, W. J., ... Isaacs, C. (2018). Prevalence of Homologous Recombination-Related Gene Mutations Across Multiple Cancer Types. JCO precision oncology, 2018, 10.1200/PO.17.00286. doi:10.1200/PO.17.00286

## ACKNOWLEDGEMENTS

R&D team and Product Management team, South San Francisco, CA  
R&D team, Carlsbad, CA  
R&D team, Ann Arbor, MI

